

Protocol of ER-alpha-BLA HEK293 Cell-based Assay for High-throughput Screening

DOCUMENT: ER-alpha-BLA_TOX21_SLP_Version1.0

TITLE: Protocol of ER-alpha-BLA HEK293 Cell-based Assay for High-throughput Screening

ASSAY REFERENCES:

Assay Target	Cell Lines	Species	Tissue of Origin	Assay Readout	Assay Provider	Toxicity Pathway
Estrogen receptor alpha: LBD (Recombinant)	HEK293	Human	Embryonic kidney	Beta lactamase reporter	Invitrogen	NR signaling

QUALITY CONTROL PRECAUTIONS:

1. -Cell culture is maintained by passaging twice a week and should not reach more than 90% confluence
2. -Culture medium should be replaced with Assay medium overnight prior to the assay
3. -The assay should be performed in black-clear bottom 1536 well plates, so the bottom of the plates should not be touched

MATERIALS and INSTRUMENTS:

Supplies/Medium/Reagent	Manufacturer	Vender/Catalog Number
Phenol red-free DMEM	Invitrogen	Invitrogen/21063
DMEM	Invitrogen	Invitrogen/11965
Dialyzed FBS	Invitrogen	Invitrogen/26400
Charcoal stripped FBS	Invitrogen	Invitrogen/12676
NEAA	Invitrogen	Invitrogen/11140
Sodium pyruvate	Invitrogen	Invitrogen/11360
Penn-strep	Invitrogen	Invitrogen/15140
Hygromycin B	Invitrogen	Invitrogen/10687
Zeocin	Invitrogen	Invitrogen/R25001
Black-clear bottom 1536 well plates	Greiner	Greiner/789092F
LiveBLAzer B/G FRET substrate	Invitrogen	Invitrogen/K1028
Recovery Cell Culture Freezing	Invitrogen	Invitrogen/12648

Medium		
0.05% Trypsin-EDTA	Invitrogen	Invitrogen/25300
Envision Plate Reader	Perkin Elmer	Perkin Elmer
BioRAPTR FRD dispenser	Beckman Coulter	Beckman Coulter
Multidrop COMBI	Thermo Electron Corporation	Thermo Electron Corporation
Beta-Estradiol (Agonist control compound)	Sigma	Sigma/E8875

PROCEDURE:

1. Cell handling:

1.1. Media Required:

Component	Growth Medium	Assay Medium	Thaw Medium	Freezing Medium
Phenol red-free DMEM	-	98%	-	-
DMEM	90%	-	90%	-
Dialyzed FBS	10%	-	10%	-
Charcoal stripped FBS	-	2%	-	-
NEAA	0.1mM	0.1mM	0.1mM	-
Sodium pyruvate	1mM	1mM	1mM	-
Penn-strep	100U/ml-100ug/ml	100U/ml-100ug/ml	100U/ml-100ug/ml	-
Hygromycin B	80ug/ml	-	-	-
Zeocin	100ug/ml	-	-	-
Recovery Cell Culture Freezing Medium	-	-	-	100%

1.2. Thawing method

1.2.1 1ml frozen cells of ERalpha-bla were taken in pre-warmed 10ml of thaw medium for centrifuging.

1.2.2 2-3ml of the thaw medium is taken to resuspend the pellet

1.2.3 The cells were seeded in T-75 flask at 2 million cells

1.3. Propagation method

1.3.1 The cells are detached using 0.05% Trypsin

1.3.2 The cells are further passaged at a density of 4-5 million cells per T-225 flask

2. Assay Protocol

- 2.1 Rinse the cells with DPBS and detach them by using 0.05% Trypsin and centrifuge
- 2.2 Resuspend the pellet with assay buffer
- 2.3 Plate the cells in black-clear bottom 1536 well plate at 5000/well/6uL through 8 tip of a plate dispenser (Multi drop)
- 2.4 Incubate at 37C for 5hrs
- 2.5 Transfer 23nL of the compounds from the library collection and positive control through Pintool
- 2.6 Incubate at 37C for 18hrs
- 2.7 Add 1uL of CCF4 dye using a single tip of a plate dispenser (Bioraptr)
- 2.8 Incubate at room temperature for 2hrs
- 2.9 Read the fluorescence intensity through Envision plate reader

3. Assay Performance

ERα-bla (Beta-Estradiol; Agonist control)	Online Validation Agonist (Mean \pm SD)
EC50	0.40 \pm 0.07 nM (n = 27)
S/B	3.68 \pm 0.19
CV (%) [*]	10.04 \pm 1.02 (n = 18)
Z'	0.73 \pm 0.05

*CV values shown represent average of DMSO plates and low concentration plates only.