

# Protocol of PGC/ERR Hek293T Cell-based Assay for High-throughput Screening

**DOCUMENT:** PGC/ERR\_TOX21\_SLP\_Version1.0

**TITLE:** Protocol of PGC/ERR Hek293T Cell-based Assay for High-throughput Screening

## ASSAY REFERENCES:

Assay Target	Cell Lines	Species	Tissue of Origin	Assay Readout	Assay Provider	Toxicity Pathway
PGC/ERR regulated genes	PGC/ERR Hek293T	Human	Embryonic kidney	Luminescence	NTP	Energy homeostasis

## QUALITY CONTROL PRECAUTIONS:

1. The cells should not be grown more than 80-85% confluence
2. Do not leave cells in Trypsin for more than 5 min at RT

## MATERIALS and INSTRUMENTS:

Supplies/Medium/Reagent	Manufacturer	Vender/Catalog Number
DMEM, high glucose	Invitrogen	11965-092
HyClone® FBS	Hyclone	SH30071-03
L-glutamine,200mM	Invitrogen	25030-081
Sodium pyruvate solution, 100 mM,	Invitrogen	11360-070
Penicillin/Streptomycin (antibiotic)	Invitrogen	15140
Puromycin	Invitrogen	A11138-03
Recovery Cell Culture Freezing Medium	Invitrogen	12648
ONE-Glo Luciferase Assay	Promega	E6051
CellTiter Fluor Cell Viability Assay	Promega	G6080
ViewLux Plate Reader	Perkin Elmer	-
Multidrop	Thermofisher	-
BioRAPTR	Beckman Coulter	-

## PROCEDURE:

1. Cell handling:
  - 1.1. Media Required:

Component	Growth Medium	Assay Medium	Thaw Medium	Freezing Medium
DMEM, high glucose	90%	90%	90%	-
HyClone® FBS	10%	10%	10%	-

Lglutamine,200mM	4mM	4mM	4mM	-
Sodium pyruvate solution, 100 mM	1mM	1mM	1mM	-
Penicillin/Streptomycin (antibiotic)	1%	1%	1%	-
Puromycin	1ug/mL	1ug/mL	-	-
Recovery Cell Culture Freezing Medium	-	-	-	100%

### 1.2. Thawing method

- 1.2.1 Place 14 mL of pre-warmed thaw medium into a 15 mL conical tube
- 1.2.2 Remove the vial of cells to be thawed from liquid nitrogen and thaw rapidly by placing at 37°C in a water bath
- 1.2.3 Mix the entire content of the vial to 14 ml of pre-warmed medium and centrifuge to remove DMSO
- 1.2.4 Discard the supernatant and reconstitute the pellet using pre-warmed thaw medium
- 1.2.5 Transfer the cells to T75 flask

### 1.3. Propagation method

- 1.3.1 Detach the cells from the flask using Trypsin-EDTA (0.05%)
- 1.3.2 The cells in growth medium are re-seeded in T225 flask

## 2. Assay Protocol

- 2.1 Harvest cells and re-suspend in assay medium and adjust the required cell density
- 2.2 Dispense 2500 cells/5  $\mu$ L/well into 1536 well tissue treated white plates using a Multidrop dispenser
- 2.3 Incubate the plates at a 37°C, 5% CO<sub>2</sub> incubator for 6 h
- 2.4 Transfer 23nL of compounds and positive control to the assay plate by a pin tool.
- 2.5 Incubate the assay plates at 37°C, 5% CO<sub>2</sub> for 17.5 h
- 2.6 Add 1  $\mu$ L of CellTiter-Fluor reagent using a BioRAPTR
- 2.7 Incubate at 37°C, 5% CO<sub>2</sub> for 17.5 h
- 2.8 Read the fluorescence intensity in ViewLux plate reader
- 2.6 Add 4uL of OneGlo and incubate the plates at room temperature for 0.5 h
- 2.7 Read the luminescence intensity in ViewLux plate reader

## 3. Assay Performance

PGC/ERR-agonist (Genistein)	Online Validation (Mean $\pm$ SD)
EC50 ( $\mu$ M)	9.16 $\pm$ 2.31
S/B	3.48 $\pm$ 0.14
CV (%)	7.22 $\pm$ 1.45
Z'	0.65 $\pm$ 0.08
N	27

PGC/ERR-antagonist (XCT790)	Online Validation (Mean $\pm$ SD)
-----------------------------	-----------------------------------

IC50 (μM)	3.97±0.29
S/B	3.01±0.12
CV (%)	7.22± 1.45
Z'	0.55± 0.07
N	27

PGC/ERR-Viability (Tetra-Br.)	Online Validation (Mean ± SD)
IC50 (μM)	N/A
S/B	7.64±0.23
CV (%)	3.96±0.85
Z'	0.90±0.03
N	27