

# Protocol of AR-BLA HEK293 Cell-based Assay for High-throughput Screening

**DOCUMENT:** AR-BLA\_TOX21\_SLP\_Version1.0  
**TITLE:** Protocol of AR-BLA HEK293 Cell-based Assay for High-throughput Screening

## ASSAY REFERENCES:

Assay Target	Cell Lines	Species	Tissue of Origin	Assay Readout	Assay Provider	Toxicity Pathway
Androgen receptor : LBD (Recombinant)	HEK293	Rat (Androgen Receptor)	Embryonic kidney	Beta lactamase reporter	Invitrogen	NR signaling

## QUALITY CONTROL PRECAUTIONS:

1. -Cell culture is maintained by passaging twice a week and should not reach more than 90% confluence
2. -The assay should be performed in black-clear bottom 1536 well plates, so the bottom of the plates should not be touched

## MATERIALS and INSTRUMENTS:

Supplies/Medium/Reagent	Manufacturer	Vender/Catalog Number
-DMEM, high glucose	-Invitrogen	-Invitrogen/11965
-Opti-MEM	-Invitrogen	-Invitrogen/11058
-Dialyzed FBS	-Invitrogen	-Invitrogen/26400
-HEPES	-Invitrogen	-Invitrogen/15630
-NEAA	-Invitrogen	-Invitrogen/11140
-Sodium pyruvate	-Invitrogen	-Invitrogen/11360
-Penicillin and Streptomycin	-Invitrogen	-Invitrogen/15140
-Hygromycin	-Invitrogen	-Invitrogen/10687
-Zeocin	-Invitrogen	-Invitrogen/R250-01
-0.05% Trypsin-EDTA	-Invitrogen	-Invitrogen/25300
-Recovery Cell Culture Freezing Medium	-Invitrogen	-Invitrogen/12648
-Black-clear bottom 1536 well plates	-Greiner	-Greiner/789092F

-LiveBLAzer B/G FRET substrate	-Invitrogen	-Invitrogen/K1028
-Multidrop COMBI	-Thermo Electron Corporation	-Thermo Electron Corporation
-BioRAPTR FRD dispenser	-Beckman Coulter	-Beckman Coulter
-Envision Plate Reader	-Perkin Elmer	-Perkin Elmer
R1881 or Methyltrienolone (Agonist control compound)	-Perkin Elmer	- Perkin Elmer/ NLP005005MG

## PROCEDURE:

### 1. Cell handling:

#### 1.1. Media Required:

Component	Growth Medium	Assay Medium	Thaw Medium	Freezing Medium
-DMEM, high glucose	-90%	-	-90%	-
-Opti-MEM	-	-90%	-	-
-Dialyzed FBS	-10%	-10%	-10%	-
-HEPES	-25mM	-	-25mM	-
-NEAA	-0.1mM	-0.1mM	-0.1mM	-
-Sodium pyruvate	-1mM	-1mM	-1mM	-
-Penicillin and Streptomycin	-100U/ml and 100ug/ml	-100U/ml and 100ug/ml	-100U/ml and 100ug/ml	-
-Hygromycin	-80ug/ml	-	-	-
-Zeocin	-80ug/ml	-	-	-
-Recovery Cell Culture Freezing Medium	-	-	-	-100%

#### 1.2. Thawing method

1.2.1 -1ml frozen cells of AR-bla were taken in pre-warmed 10ml of thaw medium for centrifuging

1.2.2 -Thaw medium is used to re-suspend the pellet

1.2.3 -Seed the cells at 2 million per T-75 flask with thaw medium

#### 1.3. Propagation method

1.3.1 -Detach the cells from the flask using 0.05% Trypsin

1.3.2 -The cells are re-seeded in T-225 flask at 3-4 million

## 2. Assay Protocol

2.1 -Spin down the cells after rinsing the cells with DPBS and trypsinizing

- 2.2 -Resuspend the pellet with assay medium
- 2.3 -Plate the cells in black-clear bottom 1536 well plate at 2000/well/6uL through 8 tip Multidrop plate dispenser
- 2.4 -Incubate at 37C for 5hrs
- 2.5 -Transfer 23nL of compounds from the library collection and positive control to the assay plates through Pintool
- 2.6 -Incubate at 37C for 16hrs
- 2.7 -Add 1uL of CCF4 (FRET Substrate) dye using a single tip plate dispenser (Bioraptr)
- 2.8 -Incubate at room temperature for 2hrs
- 2.9 -Read the fluorescence intensity through Envision plate reader

### 3. Assay Performance

<b>AR-bla (R1881; Agonist control)</b>	<b>Online Validation Agonist (Mean <math>\pm</math> SD)</b>
EC50	1.06 $\pm$ 0.10 nM (n = 27)
S/B	2.07 $\pm$ 0.18
CV (%)*	6.67 $\pm$ 2.36 (n = 18)
Z'	0.34 $\pm$ 0.13

\*CV values shown represent average of DMSO plates and low concentration plates only